







Abletech Orthopedics Corporation is easily accessible orthopedic support manufacturer and exporter, providing qualified braces to keep you moving every day for your families of all ages.

Wellcare, a house brand, is well-built with commercial package and promotional tools at affordable prices for those who have visions of this industry and those who have wisdom to distribute in their region.

Wellcare keeps you moving!

A wide range of material for daily life employed with therapeutic criteria under strict production controls plus functional lab-test accredits with ISO 13485, CE, FDA certificates which guarantee you our commitments.

Constant research, innovation and improvement on technology demanded and requested by our current working partners welt with our professional expertise in medical field which enables us being the safest and most reliable orthopedics provider globally. **Abletech Orthopedics Corporation**, without a doubt, is the perfect choice to partner.



Wellcare-supports.com


Content

NEW PRODUCT	P04	
WELLBIND SERIES	P06	
WELLKNIT SERIES	P08	
POST-OP PROTECTOR	P10	
PRESSURE CARE & RESTRAINTS	P12	
FUNCTIONAL FOOT CARE		
FITWELL SERIES		
WELLDRY SERIES		
WRAP-AROUND SERIES		
KNITTING SERIES		
KNEE BRACE	P20	
ANKLE BRACE	P26	
BACK BRACE	P32	
SHOULDER BRACE	P38	PART
ELBOW BRACE		
WRIST BRACE		
NECK BRACE		
FOOT CARE		
SIZE CHART		






NEW PRODUCT






From the medial end of the top of thigh frame, have the strap go around

thourgh the triangular buckle of upper medial calf frame. Fasten the calf strap in the same way as thigh strap.









Better strapping system for 3 point leverage offlading





42301 WellBind THUMB BRACE

SIZE S/M L/XL

- De Quervain syndrome
- Thumb ligament sprain

- · Universal design for left and right hand
- · Designed malleable aluminum stay provides stabilization
- Supportive foam and fabric provide high level comfort



41301 WellBind WRIST BRACE

SIZE S M L XL

- Carpal tunnel syndrome
- · Wrist joint sprain or overuse syndrome

- · Universal design for left and right hand
- Removable plastic stays and anatomically contoured stray enhance
- Perforated straps and foam facilitate airflow, improve compliance and comfort





31301 WellBind ELBOW BRACE WITH SILICONE

- Tennis elbow
- · Golfer's elbow

- Universal design fits left and right, one size fits all
- Removable silicone padding applies pressure relieving around painful site
- Perforated straps and foam facilitate airflow, improve compliance and comfort









WellKnit Series

For those who are allergic to latex or neoprene, Wellknit will provide superior support, proprioceptive feedback and facilitate muscle function, 3D weaving design improves breathability and flexibility. Wellknit series is ideal for symptomatic support while engaging in active mobility.



23401 WellKnit SILICONE **BACK SUPPORT**

SIZE S M L XL XXL

· Low back pain

- Removable silicone pad helps to aid firm stabilization
- WellKnit fabric fits your back snugly





23402 WellKnit BACK **STABILIZER**

SIZE S M L XL XXL

• Mild disc/ degenerative problems

- 4 plastic stays support back and help realign posture
- · WellKnit fabric fits snugly, plus side straps for better stabilization



WellKnit Series



SIZE S M L XL

Wrist instability or chronic pain

- Silicone padding for release excessive pressure around bony structures
- · Removable strap could be adjusted for additional stabilization



Palmar Side



31401 WellKnit ELBOW **BRACE**

SIZE S M L XL

• Overuse injury or ligament sprain

- Silicone padding for release excessive pressure around bony structures
- Removable strap for pressure adjustment over painful elbow joint



52401 WellKnit KNEE **BRACE**

SIZE S M L XL

INDICATIONS & FEATURES

- Perfect for chronic swelling knee or mild osteoarthritis and other mild knee condition
- Silicone padding around kneecap for protection and pressure relief
- Removable elastic upper and lower straps for stabilization



62401 WellKnit ANKLE **BRACE**

SIZE S M L XL

INDICATIONS & FEATURES

- Silicone padding relives excessive pressure around bilateral malleolus
- Removable and elastic strap fits left and right ankle, provides lateral-tomedial tension for inversion injury (also inversion tension for eversion sprain)







21101 ARM ABD ORTHOSIS

- · Following rotator cuff repairs and post operation
- Dislocation problems

- 0°, 30°, 60°, 90° immobilization
- Combined with padding and double strapping system for secure fit and comfort
- Featured arm stay for sufficient support and preventing circulatory problems





13101 MINERVA CTO

SIZE UNIVERSAL

- Post surgery or stable fractures
- Degenerative joint disease
- C3 to upper thoracic flexion injury

- Minerva type cervical thoracic orthosis
- One size fits all, adjustable with simple hex wrench and screw
- Countered frame supports occiput and chin aids in stabilization
- Easy to fit in supine position with minimal disturbance
- Removable foam padding and strapping for comfort and optimized fitting



26101 HIP ORTHOSIS

SIZE REGULAR L/R

- Following hip replacement surgery
- · After a hip dislocation or severe hip condition

- Anatomically contoured shell for better fitting
- · Adjust width and height by hex wrench and screw driver, also ensure secure application
- · Removable foam padding and strapping for comfort and optimized fitting
- Fit left and right side by switching the shaft of the angle adjustment mechanism to the other side
- · Angle adjustment mechanism: adjust flexion/extension in 20° increments, adduction/abduction in 16° increments

Post-OP Protector





Pressure Care & Restraints

For patients' safety and skin care

62801 FOOT POSITIONER

SIZE UNIVERSAL - PAIR

 Wavy and vented PU foam protect heels from skin breakdown and helps distributing pressure

• Supporting PP plate under calf muscles for better positioning

· Easy to apply with hook and loop fasteners



62802 HEEL PILLOW

SIZE UNIVERSAL - PAIR

- · Cotton filling for supporting good potision of low leg and distributing
- · Soft touch and skin-firendly liner for compatibility
- One size fits most
- · Heel opening design for easy inspection and avoiding excessive pressure







71801 WEDGE POSITIONER

SIZE UNIVERSAL

- Firm density PU foam for sleep positioning and support
- With easy-to-clean fabric lining



71802 LIMB ELEVATOR

SIZE UNIVERSAL

- Firm density PU foam for supprting arm joints or ankle joints especially for bedside care
- With easy-to-clean fabric lining



71803 JOINT POSITIONER

SIZE UNIVERSAL

- Firm density PU foam for supporting joints especially knee
- · With easy-to-clean fabric lining
- · Durable strap for secure fit







71811 LIMB HOLDER

SIZE UNIVERSAL - PAIR

- · Soft cushioning and breathable padding adds comfort
- Durable straps with quick release buckle
- No lartex

71812 HAND HOLDER

SIZE UNIVERSAL – PAIR

- Breathable hand holder with zip opening for easy inspection of fingers and oxymeter application
- Durable and light touch polyester straps

71813 CHEST RESTRAINT

- · Soft cushioning and breathable padding adds comfort
- · Durable straps with quick release buckle
- No latex

Functional Foot Care Products

Functional socks for stress relief plus gentle and efficient pain control Diabetic socks as daily partners for all day comfort Advanced Dropfoot Lifter and Arch Lifter for structural support and improvement of quality of life



61601 PLANTAR FASCIITIS SOCK

SIZE S/M L/XL

- · Calf muscle tightness and inflexibility
- Plantar fasciitis and related heel pain

- Gently stretch plantar fascia and calf muscle
- Suitable for night use: efficient pain relieving effect
- Adjust desired and comfortable pressure as needed
- Packed per piece



61602 PLANTAR FASCIITIS **COMPRESSION SOCK**

SIZE XS S M L XL - PCS

- · Calf muscle tightness and inflexibility
- Plantar fasciitis and related heel pain

- Compression sock improves circulation especially for calf mucle swelling from exercise or daily activity
- · Suitable for night use: efficient pain relieving effect
- Adjust desired and comfortable pressure as needed



61691 DIABETIC ANTIBACTERIAL SOCKS

S/M L/XL

- · Diabetic foot
- Foot conditions with wound dressings

- Bamboo charcoal and silver fiber for odor control and germproof effect
- · Extra room for toes, no excessive pressure for sensitive foot



61692 DIABETIC DAILY SOCKS

SIZE S M L

- Diabetic foot
- · Sensitive foot

- Seamless design for daily use
- · No excessive pressure for sensitive foot



62040 DROPFOOT LIFTER

SIZE UNIVERSAL

- Mild drop foot
- · Add ankle dorsiflexion lifting

- Put inside casual shoes or sneakers
- · Breathable and cushioning fabric avoid excessive friction
- Durable strap supports well



63051 ARCH LIFTER

- Supports fallen arch
- · Mild plantar factiitis

- Fits in most shoe
- · Terry cloth foam padding for lifting arch gently, elastic unapped fasteners add comfort and secure closure



FitWell Series

Innovative Design- Low Profile, Lightweight, Breathable, Fit Well.



23601 SACROFIT BELT

SIZE S M L XL

- Sacroiliac joint pain
- · SI joint hypermobility/ instability

- Stabilizing sacroiliac joint
- Help to promote injured tissue healing
- Soft and low profile air mesh plus breathable elastic band fit under clothing





23602 LUMBOFIT SUPPORT

SIZE S M L XL

- · Mild disc problems
- · Chronic back pain

- Help to stabilize low back
- Extra wide elastic straps for compressive effect and stabilization
- · Plastic stays beside avoid direct pressure on bony process and provide additional support
- Soft and low profile air mesh plus breathable elastic band fit under clothing



23603 LUMBOFIT PANEL **BRACE**

- Mild to moderate disc problems
- · Non-specific low back pain

- Back panel helps to provide stabilization and support
- Effective side straps for compressive effect, ensure better fitting
- · Low profile mesh and elastic band for breathability and comfort





Comparison

WellDry [™] Material	Traditional Neoprence	
Water-based PU	Chloroprene rubber, CR	
Cool and Ventilative	Muggy and stuffy	
Quick-drying during activity	Moist during activity	
Light and flexible	Heavy and inflexible	
Non-toxic and latex-free	Allergic and irritation	
Degradable and eco-friendly	Not degradable and pollutant	









21012 SHOULDER SUPPORT 31013 ELBOW SUPPORT

SIZE UNIVERSAL

- · Shoulder bursitis
- Shoulder tendinitis
- · Mild shoulder instability

- Ensure proper fit via closures at the arm and chest
- Help to stabilize shoulder joint in space
- · One size fits all

· Joint effusion and swelling

SIZE UNIVERSAL

Collateral ligament sprain

- · Easy to put on and take off with wraparound style
- Provide pressure relief around inflammated soft tissues

31014 ELBOW STRAP

SIZE S M L XL

- Tennis elbow (Lateral epicondylitis)
- Golfer's elbow (Medial epicondylitis)

- Silicone pad helps to provide pressure relief around inflammated soft tissues
- Adjustable hook and loop fasteners for compressive effect

Silicone Pad





SIZE S M L XL XXL

- Low back pain
- Mild sponlylosis

- Four metal stays help to provide stabilization and support
- Effective side straps for compressive effect, ensure better fitting





62039 CONTOUR ANKLE **STABLIZER**

SIZE SMLXLL/R

- Ankle instability
- Mild to moderate ankle sprain

- · Anatomically contoured stay helps to stabilize ankle joint without reducing mobility
- Extended strap adds additional support for the sprained joint





62029 ANKLE SUPPORT

- · Articular effusion and swelling
- Moderate ankle sprain and strain

- · Easy to put on and take off with wraparound style
- · Bilateral crisscross straps reinforce ankle stability and offer additional compression



41002 WRISTLET BRACE

SIZE UNIVERSAL

- · Joint effusion and swelling
- Slight wrist instability

- · Help to stabilize wrist joint
- Provide stable support and pain relief



42024 WRIST BRACE

SIZE UNIVERSAL

- · Wrist sprain and strain
- Wrist instability

- Provide comfortable compression and support to the wrist joint
- Low cut design allows for free hand movement



42025 THUMB BRACE

SIZE UNIVERSAL

- De Quervain syndrome
- Gamekeeper's thumb
- Thumb sprain and strain

- · Contoured aluminum splint immobilizes thumb movement
- Ergonomic design plus hook and loop fasteners ensures easy application and precise fit



52030 KNEE SUPPORT WITH STAYS

SIZE S M L XL

- Collateral ligament injuries
- Knee sprain and strain
- Chondromalacia

- Bilateral Spiral metal stays help to support the joint
- · Compound material design with mesh fabric provides breathability and moisture control
- · Upper and lower straps for stabilization and compression



52031 KNEE SUPPORT

SIZE UNIVERSAL

- Mild knee sprain and strain
- Mild osteoarthritis
- · Anterior knee pain

- · Easy to put on and take off with wraparound style
- Open patella design relieves pressure and enhances stabilization



52032 HINGED KNEE **SUPPORT**

- Collateral ligament injuries
- Moderate knee sprain and strain
- Meniscus injuries

- Bilateral hinges provide support and stabilization for collateral ligaments
- · Control the terminal extension angle for load relief
- Easy to put on and take off with front closure design



Wrap-around Series

- One size fits all, for left and right side
- Easy to put on and take off with wrap-around style
- Secure comfortable fitting with superior compound material
- Neoprene material provides warmth and compression, enhances circulation for healing



41601 WRAP-AROUND WRIST **SUPPORT**

SIZE UNIVERSAL

- · Joint effusion and swelling
- · Slight wrist instability

· Help to stabilize wrist joint



31008 WRAP-AROUND ELBOW **SUPPORT**

SIZE UNIVERSAL

- · Joint effusion and swelling
- Collateral ligament sprain

Provide pressure relief around inflammated soft tissues



52009 WRAP-AROUND KNEE **SUPPORT**

SIZE UNIVERSAL

- Patellar instability
- Mild knee sprain and strain
- Mild osteoarthritis

• Round patella donut offers stabilization and prevents patellar displacement



52601 WRAP-AROUND **KNEE BRACE**

SIZE UNIVERSAL

- Anterior knee pain
- Mild knee sprain

- · Short version, secure fitting, free movement
- · Spiral metal stays add support



62022 WRAP-AROUND **ANKLE SUPPORT**

SIZE UNIVERSAL

- · Articular effusion and swelling
- Mild ankle sprain

 Help to stabilize ankle joint without reducing mobility

Knitting Series

- · Knitted high-tension elastic material
- Provide warmth, breathability and moisture control
- · Support firmly and fit comfortably during movement



52019 ELASTIC KNEE **SUPPORT**

SIZE S M L XL XXL

- · Articular effusion and swelling
- Mild osteoarthritis

• Firm support, secure fitting, free movement



52020 KNEE SUPPORT

SIZE S M L XL

- · Articular effusion and swelling
- Mild knee sprain/ strain
- Mild osteoarthritis

· Spiral metal stays for additional support



52021 KNEE SUPPORT WITH SILICONE **CUSHION**

SMLXL

- Articular effusion and swelling
- Patellofemoral pain syndrome
- Mild osteoarthritis

pain relief







62023 ELASTIC ANKLE SUPPORT

SIZE S M L XL XXL

- Articular effusion and swelling
- Mild ankle sprain

- Silicone pad for pain relief and massage effect (62024)
- · Size choice for 62024: S M L XL





31010 ELASTIC ELBOW 42004 PALM BRACE **SUPPORT**

SIZE S M L XL

 Joint effusion and swelling

 Silicone pad for pain relief and massage effect (31011)



31011 with silicone



SIZE S M L XL

- Joint effusion and swelling
- · Wrist ligament sprain

 Silicone insert for pain relief and massage effect







SIZE SMLXLL/R

- Knee instability due to ACL/ PCL problem
- To restrict the range of motion after reconstructive surgery
- Recommended for low to medium impact sports

- · Adjustable ROM hinge can be locked at -10°, 0°,10°,20°,30°
- ROM control in 10° increments: Flexion -10° to 120°, Extension -10° to 30°
- Simply pull, rotate and quick-lock flexion/ extension settings
- Knee-contoured condyle pads for superior comfort
- Rigid frame and six strong straps provide desirable stability



Cousioning Padding and Updated Buckles for Comfortable and Easy Application



TRACKING STABILIZER

SIZE S M L XL

- Patellofemoral pain syndrome with maltracking or patella dislocation
- · Patellar tendonitis
- Anterior knee pain

- · Detachable design for right and left
- Patellar buttress and straps help stabilize the patella in proper position
- Ergonomical biaxial hinge provide stability and support





52603 KNEE WRAP **STABILIZER**

- · Collateral ligament and meniscus iniuries
- Moderate knee sprain and strain
- Mild to moderate osteoarthritic knee

- Open wrap design for secure and individual fit
- Double upper straps enhance support
- Light patellar buttress and wide opening for comfortable touch and prevent excessive pressure around an terior knee
- · Ergonomical biaxial hinge provide stability and support







52004 FUNCTIONAL KNEE BRACE

SIZE SMLXLL/R

- ACL / PCL / MCL / LCL and combined instabilities
- Protection and restriction of ROM after reconstructive surgery
- Recommended for low to medium impact sports

- Adjustable ROM hinge can be locked at 0°,15°,30°,60°,90°
- Cushioning and anti-slip contoured knee condyle pads hold joint in place
- Rigid frame and six strong straps provide desirable stability
- Ergonomically contoured frame padding adds comfort



52028 KNEE LIGAMENT DEFENDER

SIZE S M L XL L/R

- ACL/ PCL/ MCL/ LCL and combined instabilities
- Protection and restriction of ROM after reconstructive surgery
- Recommended for low to medium impact sports

- ROM control limitation: Flexion at 45°,60°,75°,90°;
- Extension at 0°,10°,20°,30°,40°
- Polycentric hinge with plug-in flexion/ extension stops
- High-strength streamline aluminum frame provides superior
- Soft condyle pads and padded straps for desirable comfort

















52001 **POST-OP KNEE BRACE**

SIZE UNIVERSAL

• Immobilization or ROM control for ligament and fracture-related surgeries

- Low-profile ROM hinge can be locked at 0° to 90°
- Simply rotate and quick-lock flexion/ extension settings
- ROM control in 10° increments: Flexion 0° to 120°, Extension 0° to 90°
- Telescoping design for tool-free length adjustment from 19.5" - 24.5"



52002 **POST-OP KNEE BRACE**

SIZE UNIVERSAL

• Immobilization or ROM control for ligament and fracture-related surgeries

- Low-profile ROM hinge can be locked at 0° to 90°
- ROM control in 10° increments: Flexion 0° to 120°, Extension 0° to 90°
- Simply rotate and quick-lock flexion/ extension settings
- Full wrap version with comfortable and skin-friendly foam padding material

52003 **POST-OP KNEE BRACE**

- Post-operative immobilization
- Rang of motion control for stabilization
- Following collateral/ cruciate ligaments surgery

- Tool-free length-adjustable collateral bars prevent valgus and varus deformity
- Newly designed adjustable ROM hinges with super wide range setting
- Extension limitation from 0°-120°. Flexion limitation from 0°-120°. Immobilization fixation from 0°-90°. Adjustable in every 15° increments
- · Simply rotate and press lock setting
- Easy to put on and take off by sidebuckles
- Length adjustable from 19" to 24"









52016 KNEE IMMOBILIZER

SIZE S M L XL

- Post-operative immobilization
- ACL / PCL / MCL / LCL injuries

- Three panel design for individualized
- Provide firm immobilization and complete support
- Breathable foam and skin-friendly material
- Available in lengths 16", 18", 20", 22", 24"



52007 NEOPRENE KNEE **SLEEVE**

SIZE S M L XL

- · Articular effusion and swelling
- Knee sprain and strain
- Mild osteoarthritis

- Neoprene material provides warmth and compression, enhances circulation for healing
- Anatomical cutting for comfortable fit and protection



52008 KNEE SLEEVE WITH **OPEN PATELLA**

SIZE S M L XL

- Patellar instability
- Mild knee sprain and strain
- Mild osteoarthritis

- Neoprene material provides warmth and compression, enhances circulation for healing
- · Anatomical cutting for comfortable fit and protection
- Round patellar donut offers stabilization and prevents patellar displacement



52006 PATELLAR TRACKING **ORTHOSIS**

SIZE SMLXLL/R

- · Patellofemoral pain syndrome with maltracking or patella dislocation
- · Patellar tendonitis

- Four non-elastic straps secure the orthosis against rotational forces
- Patellar buttress and straps help stabilize the patella in proper position
- Uniaxial hinge for medial/lateral support and hyperextension control



52010 PATELLAR **TENDON STRAP**

SIZE UNIVERSAL

- Patellar tendinitis (Jumpers' knee)
- Osgood-Schlatter disease

- Contoured TPE pad on patellar tendon for pain relief
- · Low-profile design for all day comfort
- Reduce strain on the patellar tendon



52034 WRAP PATELLAR **SUPPORTER**

- Patellar instability or subluxation
- Patellar lateral tilt
- Chondromalacia

- · Easy to put on and take off with front closure design
- · Additional elastic straps for compressive effect



52012 HINGED KNEE **SUPPORT**

SIZE S M L XL

- Collateral ligament injuries
- Moderate knee sprain and strain
- Patellar subluxation/ tilting
- Meniscus injuries

- Provide support and stabilization for collateral ligaments
- · Control the terminal extension angle for load relief
- · Round patellar donut relieves pressure on patella



52013 WRAP-AROUND HINGED **KNEE SUPPORT**

SIZE S M L XL

- Collateral ligament injuries
- Moderate knee sprain and strain
- Meniscus injuries

- · Provide support and stabilization for collateral ligaments
- Control the terminal extension angle for load relief
- Easy to put on and take off with front closure design



52014 POLYCENTRIC HINGED **KNEE SUPPORT**

S M L XL

- Collateral ligament injuries
- Cruciate ligament sprain
- Post meniscus surgery

- 13" length, also available for 16"(52015)
- · Controls flexion and extension at 0°, 15°, 30°, 60°, 90°
- Provide support and stabilization for collateral ligaments





52035 KNEE ROM **STABILIZER**

SIZE S M L XL

- · Mild to moderate osteoarthritis
- Mild instability

- Easy to put on and take off with lateral pulls and front open design
- · Both side stays provide stability and pressure relief
- · Breathable mesh lining and soft material provide support and comfort



52017 STABILIZED KNEE **SUPPORT**

SIZE S M L XL

- Mild collateral ligament injuries
- Knee sprain and strain
- Chondromalacia
- · Mild osteoarthritis

- Provide collateral support by spiral metal stays
- Hook and loop closures allow adjustment for ideal compression
- Round patellar donut offers stabilization and prevents patellar displacement




52018 KNEE SUPPORT WITH **MASSIVE DOTS**

SMLXL

- Patellar instability or subluxation
- · Patellar lateral tilt
- Chondromalacia

- Provide collateral support by spiral metal stavs
- Special-designed silicone insert offers pain relief and massage effect
- Comfortable wearing with mesh fabric material on the back



62034 GoFit WALKER

SIZE S M L XL

- Stable foot or ankle fracture
- Severe ankle sprain
- Post surgery application

- Bilateral air bags separately control the compression level as needed
- New air pump allows for efficient inflation and deflation
- Contoured boot shell for superior stabilization as 62011 Power Walking
- · Cushioning and lightweight rocker bottom assists natural and steady gait
- Breathable foam, skin-friendly and fast-wicking lining, soft velour fabric for maximum comfort











Designed pump protects deflation valve from misleading touch



62011 POWER WALKING BOOT

- Stable foot or ankle fracture
- Post surgery application
- Severe ankle sprain

- · Contoured boot shell provides superior stabilization
- Easy to adjust pressure with imbedded air pump system
- · Cushioning and lightweight rocker bottom assists natural and steady gait
- Liner with pad kit plus strapping ensures comfortable fit





62002 WALKING BOOT

SIZE XS S M L XL

- Stable foot or ankle fracture
- Post surgery application
- Severe ankle sprain

- Supportive but slightly malleable aluminum uprights for stabilization and fit
- · Low profile rocker bottom design assists natural and steady gait
- · Liner with pad kit plus strapping ensures comfortable fit





62004 AIR WALKING BOOT

SIZE XS S M L XL

- Stable foot or ankle fracture
- Post surgery application
- Severe ankle sprain

- Supportive but slightly malleable aluminum uprights for stabilization and fit
- Easy to adjust pressure with imbedded air pump system
- Low profile rocker bottom design assists natural and steady gait
- Liner with pad kit plus strapping ensures comfortable fit





62006 ROM WALKING BOOT

SIZE XS S M L XL

- Stable foot or ankle fracture
- Post surgery application
- Severe ankle sprain

- Control ankle motion from 45° plantarflexion to 45° dorsiflexion in 7.5° increments
- Supportive but slightly malleable aluminum uprights for stabilization
- Low profile rocker bottom design assists natural and steady gait
- Liner with pad kit plus strapping ensures comfortable fit





ANKLE SUPPORTS

62008 ROM AIR WALKING BOOT

SIZE XS S M L XL

- Stable foot or ankle fracture
- Post surgery application
- Severe ankle sprain

- Control ankle motion from 45° plantarflexion to 45° dorsiflexion in 7.5° increments
- · Supportive but slightly malleable aluminum uprights for stabilization and fit
- · Low profile rocker bottom design assists natural and steady gait
- Easy to adjust pressure with imbedded air pump system
- Liner with pad kit plus strapping ensures comfortable fit





62030 SUPER AIR WALKING BOOT

SIZE S M L XL

- Stable foot or ankle fracture
- Post surgery application
- Severe ankle sprain

- · Light weight walking boot with supportive side
- Easy to adjust pressure with imbedded air pump
- Rocker bottom design assists natural and steady
- Liner with pad kit plus strapping ensures comfortable fit





62032 SUPER WALKING BOOT

- Stable foot or ankle fracture
- · Post surgery application
- · Severe ankle sprain

- · Light weight walking boot with supportive side
- · Rocker bottom design assists natural and steady
- Liner with pad kit plus strapping ensures comfortable fit









62015 NIGHT SPLINT

SIZE S M L

- Plantar fasciitis
- · Achilles tendonitis

- Soft foam liner and padded straps provide excellent comfort
- Optional wedge for additional stretch if needed
- · Anti-slip outsole ensures steady steps and standing
- Suitable for night use: efficient pain relieving effect



62016 SOFT NIGHT SPLINT

SIZE S M L

- Plantar fasciitis
- · Achilles tendonitis

- · Skin-friendly fabric and washable, soft foam padding system
- Soft padded shell relieves pressure and provides protection
- Optional wedge for additional stretch if needed
- Suitable for night use: efficient pain relieving effect



62017 DORSAL NIGHT SPLINT

SIZE S/M L/XL

- Plantar fasciitis
- Achilles tendonitis

- Suitable for night use: efficient pain relieving effect
- Low profile and lightweight design enhances compliance
- · Plastic shell allows a gentle stretch of the plantar fascia
- · Hook and loop closures for adjustable



62012 AFO SPLINT

SIZE S M L XL L / R

• Drop foot due to hemiplegia, peroneal nerve palsy or orthopedic conditions

- Pre-fabricated polypropylene (PP) ankle-foot orthosis
- · Provide a dorsiflexion assist while walking
- · Low profile and easy to trim with scissors



62013 AFO LEAF SPRING

SIZE S M L XL L / R

• Drop foot due to hemiplegia, peroneal nerve palsy or orthopedic conditions, better for flaccid type dropfoot

- · Pre-fabricated polypropylene (PP) ankle-foot orthosis
- · Provide a dorsiflexion assist while walking
- · Lightweight and comfortable open heel design
- · Low profile and easy to trim with scissors

ANKLE SUPPORTS



62019 ANKLE SUPPORT

SIZE S M L XL

- · Articular effusion and swelling
- Mild ankle sprain

- Neoprene material provides warmth and compression, enhances circulation for healing
- · Anatomical cutting for comfortable fit and protection
- · Open-heel design for precise fit and easy application



62020 ANKLE BRACE WITH STRAP

SIZE S M L XL

- Ankle instability
- Mild to moderate ankle sprain

- Neoprene material provides warmth and compression, enhances circulation for healing
- Spiral metal stays restrict lateral movement and provide stabilization
- Elastic strap reinforces the ankle against inversion torsion
- Open-heel design for precise fit and easy application



62021 ANKLE BRACE WITH STRAP

SIZE S M L XL

- Ankle instability
- Mild to moderate ankle sprain

- Neoprene material provides warmth and compression, enhances circulation for healing
- Elastic strap reinforces the ankle against inversion torsion
- Open-heel design for precise fit and easy application



62018 FOOT LIFT

SIZE S M L

· Support drop foot raising during ambulation

- Powerful elastic strap can assist and stimulate ankle raising movement
- · Impeding normal ankle movement during ambulation
- · Low profile and breathable fabric design enhances compliance
- · Easy to put on and take off via quick release



62025 AIR STIRRUP BRACE

SIZE UNIVERSAL

- · Articular effusion and swelling
- Moderate ankle sprain with instability
- Chronic ankle instability


- · Anatomically-contoured shells support subtalar ioint
- Easy-to-adjust heel strap for precise fit and comfort
- Inflatable air cushion offers comfort and compression to reduce edema



62026 AIR/GEL STIRRUP **BRACE**

- Articular effusion and swelling
- · Acute ankle sprain with instability
- · Chronic ankle instability

- Anatomically-contoured shells support subtalar joint
- · Bladder combines air and gel for securely stabilizing the ankle joint
- · Gel can be placed in the freezer for cold therapy to swollen ankles



24001 JEWETT HYPEREXTENSION BRACE

SIZE S M L XL

- Postural kyphosis
- Stable spinal compression fractures (T7-L2) without osteoporosis

- Easy to adjust the width and height of frame
- · Lightweight aluminum alloy frame with soft padding
- Three-point pressure control provides stabilization and helps to correct kyphotic
- Flexible pelvic band and floating sternal plate for optimal fit and comfort



24002 THORACOLUMBAR SPINE SUPPORT

- Following lumbar laminectomy, discectomy or spinal fusion
- Kyphosis with chronic back pain
- Immobilization for below T7 thoracolumbar fractures

- Paraspinal uprights and chairback design maintain the spine in neutral position
- · High-tension elastic straps produce intracavitary pressure and offer stability
- Bilateral shoulder pads provide optimal comfort









23020 PULLEY SPINAL BRACE

SIZE S M L XL XXL

- Mild disc problems
- Chronic back pain
- Stable bony injury

- Enhance core muscle function with adjustable pulley system
- Keep abdominal pressure in stable for support and improve tissue healing after injuries around low back
- Extra firm and breathable foam for superior wearing comfort and support



23021 LSO PULLEY SPINAL BRACE

SIZE S M L XL XXL

- · Mild to moderate disc problems
- Chronic back pain
- Following spinal surgery

- Back panel adds stabilization
- · Enhance core muscle function with adjustable pulley system
- Keep abdominal pressure in stable for support and improve tissue healing after injuries around low back



Air Pump System



23016 AIRFOAM LUMBAR **SUPPORT**

SIZE S M L XL XXL

- Mild spondylolysis
- Mild disc problems
- Degenerative condition

- · Light and thin neoprene material plus breathable air mesh for optimal fit
- Easy to adjust pressure with imbedded air pump system
- Cross straps add stabilization around low back
- Low profile style for all day comfort

23003 ABDOMINAL BINDER

SIZE S M L XL XXL

· Postoperative fixation of abdomens

- Provides fixation to both upper and lower abdomen
- Breathable foam and skin-friendly lining for wearing comfort
- High elastic material adds firm compressive support



23004 ELASTIC LUMBAR **SUPPORT**

SIZE S M L XL XXL

- Low back pain
- · Low back sprain and strain
- Degeneration of the spine

- High elastic material adds firm compressive support
- Flexible stays plus effective side-pulls provide firm support and secure fit





23005 LUMBAR SUPPORT WITH STAYS

SIZE S M L XL XXL

- Low back pain
- Low back sprain and strain
- · Degeneration of the spine

- · Made of high quality vented elastic band for comfort and breathability
- Effective side-pulls provide firm support and secure fit
- Metal stays provide support and stabilization







23006 THERMO LUMBAR SUPPORT

Front View

SIZE S M L XL XXL

VIDICATIONIC

- Low back pain
- Low back sprain and strain
- Degeneration of the spine

FATURES

- Neoprene material provides warmth and compression, enhances circulation for healing
- Metal struts and bilateral elastic side-pulls offer optimal stability
- Pliable stays gradually assist the lumbar spine to the normal curvature



23007 LUMBOSACRAL SUPPORT

SIZE S M L XL XXL

INDICATIONS

- Back sprain and strain
- Disc problems
- Mild spondylolisthesis

FEATURES

- Foam padding material for support and comfort
- Six metal stays restrict over lordotic posture and add stabilization
- Effective side-pulls provide firm support and secure fit



23009 LUMBAR STABILIZER

SIZE S M L XL

INIDIO ATIONIC

- Back sprain and strain
- Disc problems
- Mild spondylolisthesis

FEATURES

- Additional polymer pad for pain relief and massage effect
- Air mesh material for breathability and comfort
- Keep abdominal pressure in stable for stabilization



Polymer Pad

23010 MOLDED BODY JACKET

SIZE S M L XL

- Following lumbar laminectomy or discectomy
- Below T11 thoracolumbar compression fractures



- · Removable anatomically contoured rigid anterior and posterior panel for stabilization
- · Effectively restrict movement of thoracolumbar spine
- · Hook and loop closures for secure fit and support





23011 MATERNITY SUPPORT BELT

SIZE S M L XL

· Low back pain and discomfort during pregnancy

- Lift the weight of belly and provide stabilization to low back
- Bilateral stretchable side-pulls for greater support



26022 COUNTER HIP BRACE

SIZE S M L XL L / R

- Hip sprain and strain
- Reminder for following hip orthosis removal

- Side pocket with cushioning pad helps to support and protect hip joint
- Highly breathable and skin-friendly material for optimal fit and comfort
- · Easy to put on and take off with wrap around





SHOULDER BRACE



21014 PILLOW ARM ABD **IMMOBILIZER**

SIZE S M L XL

- Following shoulder surgery
- Shoulder dislocation or subluxation
- · Shoulder hand syndrome

- ABD pillow keeps shoulder joint in place
- Hand grip for better support, shoulder strap for pressure



21018 METAL ARM ABD **IMMOBILIZER**

SIZE S M L XL

- Following shoulder surgery
- · Shoulder dislocation or subluxation
- Shoulder hand syndrome

- Adjustable metal stay keeps shoulder joint in place
- Hand grip for better support, shoulder strap for pressure



21001 SHOULDER ABDUCTION **IMMOBILIZER**

SIZE S M L XL

- · Following shoulder surgery
- Shoulder dislocation or subluxation

- Immobilize shoulder joint in 15° abduction
- · Foam padded shoulder straps for optimal comfort
- · Contoured pillow adds stability and support
- · Additional exercise ball facilitates hand function and enhances circulation



21002 SHOULDER **IMMOBILIZER AB45**

SIZE S M L XL L / R

- · Following shoulder surgery
- · Shoulder dislocation or subluxation

- Immobilize shoulder joint in 45° abduction
- Foam padded shoulder straps for optimal comfort
- · Contoured pillow adds stability and support
- · Additional exercise ball facilitates hand function and enhances circulation





21003 MESH SLING/ 21004 TC SLING

S M L XL

- · Protect and stabilize injured arm and shoulder
- Support fractured arm

- · Thumb loop and shoulder strap provide stability
- Breathable mesh lining for optimal compliance (21003)
- Made of durable, light and breathable fabric (21004)



21013 HUMERAL FRACTURE **BRACE**

SIZE S M L XL

Mid to distal third closed humeral shaft fractures

- · Ideal application to assist early and proper healing
- Adjustable straps offer desired compression to soft tissues
- Allow full range of motion at shoulder and elbow joints
- Fits left and right sides
- With two cotton stockinettes and an arm sling



21005 UNIVERSAL SWATHE **SLING**

SIZE UNIVERSAL

- Immobilization after surgical repair and humeral lesion
- Support to upper extremity paralysis and hemiplegia

- · Arm cradle style for comfortable and secure fit
- · Chest swathe enhances the stabilization of shoulder joint



21006 SHOULDER SUPPORT

- Shoulder instability
- Shoulder tendinitis or overuse injury
- · Chronic shoulder pain following hemiparesis

- Protect and support shoulder joint
- Easy to adjust with hook and loop fasteners

SHOULDER BRACE



21007 POSTURE SUPPORT **BRACE**

SIZE S M L XL

• Clavicle fractures and postural correction

- · Made of high quality vented elastic band for comfort and
- · Crisscross design adds additional trunk support



21008 POSTURE CORRECTIVE **BRACE**

SMLXL

• Clavicle fractures and postural correction

- · Made of high quality vented elastic band for comfort and breathability
- Easy to put on and take off with front closure design



21009 UNIVERSAL CLAVICLE **SUPPORT**

SIZE <u>UNIVERSAL</u>

- Clavicle fractures
- Round shoulder

- Foam padded shoulder straps with velvet lining for cushioning and stabilization
- · Easy to adjust with hook and loop fasteners



21010 CLAVICLE SUPPORT

- · Clavicle fractures
- · Round shoulder

- · Foam padded shoulder straps with velvet lining for cushioning and stabilization
- · Easy to adjust with hook and loop fasteners





31002 POST-OP ELBOW **BRACE**

SIZE UNIVERSAL L / R

- Following elbow surgery
- Stable fractures

- Control elbow flexion from 0° to 120° and extension from 0° to 90° in 10° increments
- Foam padded shoulder straps with velvet lining for cushioning and stabilization
- Adjust with simple rotate and press-lock settings
- Easy to adjust size with telescoping length adjustment design
- Specifically for left and right elbow
- Provide single brace without strap as customer's order:









31003 ELBOW PROTECTOR

SIZE S M L XL

- Tennis elbow/ golfer's elbow
- Overuse injury or ligament sprain

- Easy to adjust with hook and loop fastener
- Plastic shell adds compression and stabilization
- Provide pressure relief around inflammated soft tissues



31004 ELBOW STRAP WITH PAD

- Tennis elbow/ golfer's elbow
- Overuse injury or ligament sprain

- Imbedded padding provide pressure relief around inflammated soft tissues
- · Neoprene material provides warmth and compression, enhances circulation for healing





31005 ELBOW STRAP

SIZE S M L XL

- Tennis elbow/ golfer's elbow
- Overuse injury or ligament sprain

- Neoprene material provides warmth and compression, enhances circulation for healing
- Provide pressure relief around inflammated soft tissues



31006 AIRFORM ELBOW STRAP

SIZE UNIVERSAL

- Tennis elbow/ golfer's elbow
- Overuse injury or ligament sprain


- Inflated air cell adds pain relief and compression
- Provide pressure relief around inflammated soft tissues



31009 NEOPRENE ELBOW **SLEEVE**

- Joint effusion and swelling
- Overuse injury or ligament sprain

- Neoprene material provides warmth and compression, enhances circulation for healing
- · Air mesh fabric adds breathability and comfort







42006 FUNCTIONAL **POSITION SPLINT**

SIZE XS S M L L / R

- · Following wrist and hand surgeries
- For all cases of wrist hand dysfunction
- Neurological disabilities

- Prefabricated splint keeps wrist and hand in functional position
- Made of breathable foam and skinfriendly lycra material
- Easy to adjust with hook and loop fasteners
- · Designed specifically for left and right hand



42026 UNIVERSAL **WRIST SPLINT**

SIZE S/M L/XL

- · Carpal tunnel syndrome
- Following wrist and hand surgeries
- · Rheumatoid arthritis

- Universal design with palmar splint fits both left and right hand
- Support wrist joint in neutral position
- Breathable material with soft foam padding for all day comfort
- · Easy to fit and adjust with hook and loop closures



42027 CTS WRIST ORTHOSIS

SIZE XS S M L

- · Carpal tunnel syndrome
- Following wrist and hand surgeries
- · Immobilize wrist and fingers
- Overuse injury

- · Immobilization from forearm to distal finger joints
- · Contoured and malleable palmar splint provide optimal fit and stabilization
- Light weight low profile style for comfort and compliance
- Universal design fits both left and right hand







WRIST BRACE



42002 THUMB SPLINT

SIZE SMLL/R

- De Quervain syndrome
- · Gamekeeper's thumb
- Stabilization of 1st MP joint

- Perforated shell improves breathability
- Ergonomically contoured plastic shell for stabilization
- · Designed specifically for left and right thumb



- De Quervain syndrome
- Gamekeeper's thumb
- Thumb ligament sprain

Palmar View

- Anatomically contoured thumb splint restricts thumb movement
- · Neoprene material provides warmth and compression, enhances circulation for healing
- · Designed specifically for left and right thumb



42003 WRAP-AROUND WRIST BRACE

SIZE UNIVERSAL

- Wrist instability
- Overuse injury or wrist sprain

- · Wrap-around and universal style fits both left and right
- · Low cut design allows for free hand movement



42001 UNIVERSAL WRIST BRACE

SIZE UNIVERSAL

- · Carpal tunnel syndrome
- Overuse injury or wrist sprain

- · Universal design with removable palmar stay, fits both left and right hand
- Help to keep wrist joint in neutral, relieve pain and stress







42007 WRIST/FOREARM SPLINT

SIZE SMLXLL/R

INDICATIONS

- Carpal tunnel syndrome
- Overuse injury or wrist sprain
- Following wrist surgery

FFATURES

- Perforated suede material adds breathability
- · Longer palmar stay enhances stability



42010 WRIST SPLINT

SIZE SMLXLL/R

INDICATIONS

- Overuse injury or wrist sprain
- Carpal tunnel syndrome

FFATURES

- Breathable foam and air mesh material adds breathability and comfort
- Palmar and dorsal stays provide support and stabilization







42012 WRIST SPLINT WITH THUMB

SIZE SMLXLL/R

INDICATIONS

- · Carpal tunnel syndrome
- Overuse injury or wrist sprain
- De Quervain syndrome

FFATURES

- Palmar and thumb stays help to stabilize thumb and wrist
- Made of breathable foam and lycra material for optimal fit and comfort



42008 WRIST BRACE

SIZE SMLXLL/R

INDICATIONS

- Carpal tunnel syndrome
- Overuse injury or wrist sprain

FEATURES

- High elastic band plus a strap support and stabilize wrist joint
- Palmar stay restricts wrist movement



42009 WRIST BRACE

SIZE SMLXLL/R

INDICATIONS

- Carpal tunnel syndrome
- Overuse injury or wrist sprain

FEATURES

- High elastic band stabilizes wrist joint
- Palmar stay restricts wrist movement
- Easy to adjust with hook and loop closures





42011 WRIST BRACE

SIZE SMLXLL/R


INDICATIONS

- Carpal tunnel syndrome
- Overuse injury or wrist sprain

FEATURES

- · Palmar stay helps to stabilize wrist joint
- Made of breathable foam and lycra material for optimal fit and comfort
- Low cut design allows for free hand movement









12001/12002 SOFT COLLAR

SIZE XS S M L XL

- · Immobilization for cervical spine
- Neck sprain or strain
- Cervical spondylosis

- · Medium density foam collar with chin support
- Limit excessive neck movement
- · Easy to adjust with hook and loop fasteners
- 12002: firm density version



12003 TRACHEOTOMY COLLAR

SIZE S M L XL

- Following spinal surgery
- Stable cervical fracture

- · Lightweight and firm foam support plus plastic splint reinforcement
- With trachea opening
- Anatomical design fits chin and neck, supports cervical spine in neutral position



12004 RIGID SPLINT COLLAR

SIZE S M L XL

- Following spinal surgery
- Neck sprain or strain
- · Cervical spondylosis


- · Help to keep cervical spine in neutral
- · Perforated splint with foam padding adds comfort and
- · Easy to adjust with hook and loop fasteners



12005 RIGID SPLINT COLLAR WITH CHIN SUPPORT

- Following spinal surgery
- Cervical spondylosis

- Help to keep cervical spine in neutral
- · Perforated splint with foam padding adds comfort and compliance
- · With chin support







63001 FULL LENGTH INSOLE

SIZE XS S M L XL

• Ideal for heel pad pain and metatarsalgia. The dual-density silicone insole is capable of absorbing ground reaction force during different phases of the gait cycle. It also helps alleviate pain and provide cushion.



63002 3/4 LENGTH INSOLE

SIZE S M L

· Lower-density blue silicone enhances the shock absorption of heel-strike. With a builtin longitudinal arch and heel elevation, the orthosis is capable of absorbing up to 40% of shock impact at heel-strike. Recommended for plantar fasciitis, heel spurs and heel pain.



63003 SILICONE **INSOLE WITH** METATARSAL PAD

SIZE S M L

• Ideal for relieving forefoot pressure of high-heel user. It has a builtin soft blue dome, which helps to keep metatarsal head in natural alignment and avoids callus builtup.



63004 HEEL CUP **CENTRAL SPUR**

SIZE S M L XL

• The dual density extra soft silicone heel cup is ideal for heel pad pain or spurs formation in the middle area. It helps relieve pressure and provide additional heel stability.



63005 HEEL CUP LATERAL SPUR

SIZE S M L XL

• The dual density silicone orthosis provides heel support and additional relief to sensitive heel pad pain or heel spurs in the lateral area. Heel cushions also reduce the strain on Achilles tendon.



63006 HEEL CUSHION WITH REMOVABLE **PAD**

SIZE S M L

 It offers protection and relief to heel spur formation. The removable pad design distributes body weight away from the affected area by unloading painful spur. It also reduces strain on ligaments and Achilles tendon.



63007 HEEL CUP

SIZE S M L XL

• The silicone inserts help to protect the calcaneus from shock impact or abrasion. Heel cup increases the stability of ankle joint. Recommended for prevention of painful heel and plantar fasciitis.



63008 "U" SHAPED **HEEL CUSHION**

SIZE S M L

- It serves as a support for securing stability to the calcaneus.
- The silicone material also offers cushion function to absorb shock impact and redistribute pressure away from painful area.



64001 **BUNION CORRECTOR**

SIZE SMLL/R

• Ideal for hallux valgus. The aluminum stay helps to keep the big toe in normal alignment. Additional thickening foam padding with breathable fabric provides comfortable correction and immobilization.



64002 BUNION NIGHT SPLINT

SIZE SMLL/R

• Ideal for hallux valgus. The Night Splint keeps the big toe in normal position with gentle corrective torque. It helps straighten the big toe that is overlapping the second toe, which is usually related to bunion formation.



64003 TUBULAR TOE **PROTECTOR**

SIZE SML

• The product is made of 15cm length knitted fabric lined with mineral oil gel for protecting toes. It prevents all friction and shearing to the skin.



64004 TOE CAP

SIZE S/M L/XL

• The Toe Cap shelters the tip, dorsal and plantar aspects of the toes to protect the lesion located on the toes and nail. Silicone material eases friction on the corns, blisters, and calluses.





64005 TOE RING

SIZE S M L

- The Toe Ring, which is made of silicone material, protects the toes from friction and pressure affecting
- on blisters, calluses and corns.



64006 **METATARSAL PAD**

SIZE S/M L/XL

• Ideal for conditions such as fat pad atrophy, metatarsalgia and Morton's neuroma. It is able to alleviate painful friction and support to the metatarsal heads during walking.



64007 **HALLUX VALGUS SHIELD**

SIZE S/M L/XL

• The molded silicone orthosis absorb lateral shock and minimize friction on the bunion area. It brings constant relief of the affected area.



64008 5th TOE **BUNION SHIELD**

SIZE S/M L/XL

· The molded silicone orthosis is capable to prevent callus formation and help relieve pain and pressure over the bunion area.



64009/64010 **TOE SEPARATOR**

SIZE S M L

• The extra soft silicone Toe Separator isolate toes from rubbing and protect the sensitive toes. Recommended for overlapping toes and pinched toes.



64011 **BUNION SEPARATOR**

SIZE S/M L/XL

• The molded silicone orthosis absorbs lateral shock, minimizes friction on the bunion area and keeps toes from overlapping with each other. It brings constant relief of the affected site. Suitable for big toe or 5th toe bunion. Recommend S/M size for 5th toe bunion and L/ XL for big toe.







63103 FULL LENGTH **CUSHIONING** INSOLE

SIZE S/M L/XL

INDICATIONS & FEATURES

- For plantar fasciitis or heel pain syndrome
- Resilient foam and soft lining provides comfort and cushioning
- Dual density area for structural support and stress relief
- Also helps alleviate pain and provide cushion effect

63106 3/4 ARCH **SUPPORT**

SIZE S M L XL

INDICATIONS & FEATURES

- For moderate heel pain and unstable gait
- Semi-rigid arch support for structural control
- Provide advanced stabilization with gel cushioning around heel cup

63108 KIDS **ORTHOTICS**

SIZE 17 18 19 20

INDICATIONS & FEATURES

- For lifting and developing foot arch of kids
- Structural control and support aid in arch development
- Breathable fabric, supportive contoured design for comfort and compliance






KNEE BRACE

FUNCTIONAL KNEE BRACES


Size	S	М		XL
Circum. taken 6" above knee (IN)	15½–18½	18½–21	21-23½	23½–26½
Circum. taken 6" above knee (CM)	39 – 47	47–53	53-60	60–67
Circum. taken 5" below knee (IN)	11½–14½	14½–17	17–19½	19½–22½
Circum. taken 5" below knee (CM)	29–37	37–43	43–50	50–57



Size Chart Fit to: 52004, 52005+, 52028

KNEE BRACE & PATELLAR SUPPORT

Size	S	М	L	XL	
Knee Circum. (IN)	13–14	14–15	15–16	16–17	
Knee Circum. (CM)	33–35.5	35.5–38	38–40.5	40.5–43	



Size Chart Fit to: 52006, 52007, 52008, 52012, 52013, 52014, 52016, 52017, 52018, 52019, 52020, 52021, 52030, 52032, 52034, 52035, 52401, 52602, 52603

WALKER & ANKLE BRACE

WALKERS

Size	XS	S	M		XL
EU	Up to 35½	36-39	40-44	44-461/2	46½ +
US Men	Up to 4	4½-7	7½-10	10½-12½	12½ +
US Women	Up to 4½	5-81/2	9-11	11½-13½	13½ +

Size Chart Fit to: 62001, 62002, 62003, 62004, 62005, 62006, 62007, 62008, 62011, 60030, 60031, 60032, 60033, 62034

ANKLE BRACES & SUPPORTS

Size	S	М	L	XL	
Ankle Circum. (IN)	6–8	8–10	10–12	12–14	
Ankle Circum. (CM)	15–20	20–25	25–30	30–35	



Size Chart Fit to: 62019, 62020, 62021, 62023, 62024, 62039, 62401

FOOT LIFT

Size	S	М	L	
Ankle Circum. (IN)	7–8	8½-10	10½–13	
Ankle Circum. (CM)	18–21	22-26	27–33	



Size Chart Fit to: 62018, 62040

WALKER & ANKLE BRACE

AFO SPLINT Size 37-40 40-44 44-47 47 + US Men 5-7 7½-10 101/2-13 13 + 6-8 US Women 8-11 11½-14 14 +



Size Chart Fit to: 62012, 62013

NIGHT SPLINTS

Size	S	М	
EU	Up to 39	40-44	44 +
US Men	Up to 7	7½-10	10½ +
US Women	Up to 7½	8-101/2	11+

Size Chart Fit to: 62015, 62016

DORSAL NIGHT SPLINTS

Size	S/M	L/XL
EU	37-43	43-481/2
US Men	5-9	9½-14
US Women	6-10	10½ +

Size Chart Fit to: 62017

BACK BRACE

BACK SUPPORT & HIP BRACE

Size	S	М	L	XL	XXL
Waist Circum. (IN)	22-28	27-33	32-38	37-43	42-48
Waist Circum. (CM)	56-71	69–84	81–96	94-109	107–122

Size Chart Fit to: 23003, 23004, 23005, 23006, 23007, 23009, 23010, 23014, 23016, 23020, 23021, 23401, 23402, 23601, 23602, 23603, 24002, 24003, 24004, 26022



MATERNITY SUPPORT BELT

Size	S	М		XL	
Gestational Age (week)	3-10	11-14	16-18	20 +	

Size Chart Fit to: 23011

HIP ORTHOSIS

Size	Regular size	
Thigh circum.	50-70cm	
Waist circum.	56–105cm	

Size Chart Fit to: 26101

JEWETT HYPEREXTENSION BRACE

Size	S	М	L	XL
Pelvic Circum. (IN)	24–29	29–35	35–41	41–45
Pelvic Circum. (CM)	60–75	75–90	90–105	105–115
Trunk Length (IN)	15–18	16–19	18–20	19–22
Trunk Length (CM)	38–46	41–48	46–51	48–56

Size Chart Fit to: 24001



NECK BRACE

CERVICAL COLLAR								
Size	XS	S	М	L	XL			
Neck Circum. (IN) Neck Circum. (CM)					20–22 51–56			
Height (IN) Height (CM)	2½ 6.4	3 7.6	3½ 8.9	4 10.2	4½ 11.4			

Size Chart Fit to: 12001, 12002

TRACHEOTOMY COLLAR

Size	S	M		XL	
Neck Circum. (IN)	10"-13"	13"-16"	16"-19"	19"+	
Neck Height		~ 4¼"	~		

Size Chart Fit to: 12003

MINERVA CTO

Size	Universal	
neck circum.	36–42cm	
Size Chart Fit	to: 13101	



RIGID SPLINT COLLAR

Size	S	М	L	XL
Neck Circum. (IN)	10–11	12–14	14–16	16–18
Neck Circum. (CM)	25–28	30–36	36–41	41–46
Height (IN)	2½	3	3½	4
Height (CM)	6.4	7.6	8.9	10.2
Size Chart Fit to: 120	04, 12005			

WRIST BRACE

FUNCTIONAL POSITION SPLINT

Size	XS	S	М	L	
Measure the width of MP joint (IN)	21/4-23/4	2¾-3¼	31/4-31/2	3½-4	
Measure the width of MP joint (CM)	6–7	7–8	8-9	9–10	

Size Chart Fit to: 42006

WRIST / THUMB SUPPORT


Size	XS	S	М		XL	
Wrist Circum. (IN)	under 5	5-6	6-7	7–8	8-9	
Wrist Circum. (CM)	under 13	13-15	15 – 18	18 -20	20-23	

Size Chart Fit to: 42002, 42004, 42005, 42006, 42007, 42008, 42009, 42010, 42011, 42012, 42013, 42027, 41301, 42301, 42401

UNIVERSAL WRIST SPLINT

Size	S/M	L/XL
Wrist Circum. (IN) Wrist Circum. (CM)	5–7 13 – 18	7–9 18- 23

Size Chart Fit to: 42026




SHOULDER BRACE

SHOULDER IMMOBILIZER & SLING

Size	S	М	XL
Forearm Length (IN) Forearm Length (CM)			

Size Chart Fit to: 21001, 21002, 21003, 21004, 21014, 21018

ARM ABD ORTHOSIS



Size Chart Fit to: 21101



SHOULDER & CLAVICLE SUPPORT

Size	S	М	L	XL
Chest Circum. (IN)	24-30	30-36	36-42	42-48
Chest Circum. (CM)	61-76	76-91	91-107	107-122

Size Chart Fit to: 21006, 21007, 21008, 21010



HUMERAL FRACTURE BRACE

Biceps Circum. (IN) 8 Biceps Circum. (CM) 2		

Size Chart Fit to: 21013



ELBOW BRACE

ELBOW STRAP & SUPPORT

Size	S	М	L	XL	
Circum. of widest part of forearm (IN)	8–9	9–10	10–11	11–12	
Circum. of widest part of forearm (CM)	20–23	23–25	25–28	28–30	

Size Chart Fit to: 31003, 31004, 31005, 31009, 31010, 31011, 31014, 31401



DIABETIC CARE

DIABETIC DAILY SOCKS

Size	S	М	L
EU	36-39	39–42	42-45
US men	4-61/2	6½-9	9-111/2
US women	5-71/2	7½-10	10-121/2
Mondo	22-24.5	24.5-27	27-29.5

Size Chart Fit to: 61692

DIABETIC ANTIBACTERIAL SOCKS

Size	S/M	L/XL	
EU	37–42	42½-48	
US men	5-9	9½-14	
US women	6-10	10½-15	
Mondo	23-27	27.5-32	

Size Chart Fit to: 61691

46-48



FOOT CARE

FOOT CARE

Size	XS	S	М		XL
EU	34	36	38	40	42
US men	2½	4	5½	7	9
US women	3½	5	6½	8	10

Size Chart Fit to: 63001

Size	S	М	
EU	34-38	38-42	42-46
US men	21/2-51/2	5½-9	9-12
US women	31/2-61/2	6½10	10-13

Size Chart Fit to: 63002, 63003, 63006, 63008, 64001, 64002, 64003, 64005, 64009, 64010

PLANTAR FASCIITIS SOCKS

Calf Circum. (IN) under 16	16-21
Calf Circum. (CM) under 40	40-53

Size Chart Fit to: 61601

Size	S	М		XL
EU	35-37	37.5-40	40.5-42	42-45
US men	3½-5	5½-7	7½-9	9½-11½
US women	4½-6	6½-8	8½-10	10½-12½

Size Chart Fit to: 63004, 63005, 63007

Size	5/ I ^V I	L/ XL
EU	35-40	45
US men	3½-7	7½-11½
US women	4½-8	8½-12½

Size Chart Fit to: 64004

Calf Circum. (cm) 30-39 33-41 36-45 39-48 42-51 Ankle Circum. (cm) 18-22 21-25 26-29 30-33 33-37

Size Chart Fit to: 61602

INSOLES & ORTHOTICS

Size	S/M	L/XL	Size	S	М	L
US women	3½-7½	41–46 7½–12 8½–13 25–30	US women	6½-8½	7½-9½	43-45½ 9½-12 10½-13 27.5-30

12-14 13-15 30-32

Size Chart Fit to: 63103 Size Chart Fit to: 63106

Size	17	18	19	20
EU Kids	27½-28½	28½-30	30-311/2	31½-33½
US kids	11-12	12-13	13-1	1.5-2
CM	17.5	18.5	19.5	20.5

Size Chart Fit to: 63108

KEEPS YOU MOVING



Abletech Orthopedics Corporation

No. 8, Lane 234, Xinshu Road, Xinzhuang District, New Taipei City 242, Taiwan TEL: +886-2-2926-1765

TEL: +886-2-2926-1765 FAX: +886-2-2926-2933

E-mail: sales@wellcare-supports.com

Copyright © 2019 by Abletech Orthopedics Corporation